01 / STARTING POINT

Finding a more consistent fold

One eyelid has an established crease. The other needs a clearer starting path. This sample shows how the technician compares the two.

The goal in this case

The right eyelid is the reference. The left eyelid is the focus: find a cut size and placement that help it follow a fold closer to the right.

LOOKING SLIGHTLY DOWN



Notice the faint line above the lashes

The submitted video shows a faint line becoming easier to see when the gaze moves slightly below eye level. Compare this with the straight-ahead view.

LOW MICROCREASES



A useful similarity

Both eyelids show small existing lines very close to the upper lashes. Their low position and overall direction are broadly similar.

02 / WHAT IS DIFFERENT

Look above the lowest line

The lowest microcreases have similarities. The pattern above them is where the two eyelids differ more clearly.



What to compare in the marked areas

Compare the direction of the small lines and the spaces between them. Look at whether the skin gathers into one main fold or follows several nearby lines.

Reference eyelid: an established path

The right eyelid already gathers into a main crease. That existing fold provides the reference for this comparison.

Eyelid being reviewed: competing patterns

The left eyelid does not yet gather in the same way. The aim is a more consistent intended path rather than alternating between nearby wrinkles.

What this image does not tell us

A visible line is an observation from this case. It does not establish how long a fold will last or promise a permanent result.

03 / RESPONSE TO TRACING

A temporary fold appears

These frames are from the submitted tracing video. They show a temporary response while the eye is opened actively and the gaze moves upward.



Video frame 1



Video frame 2



Video frame 3



Video frame 4

What to notice

Look for the skin folding inward, beyond a surface line alone. In this video, the fold remains briefly while the eye is held actively open.

How this informs the next step

This temporary response gives the technician something to compare during the initial size trial. It is not evidence that the fold has already become lasting.

04 / ONE TECHNIQUE CHANGE

Use the single-point end

For this client, the recommendation is to use the single-point end rather than the two-prong fork when following the intended line.



Follow one intended path

Use controlled movements to follow the desired crease path and height. Avoid stimulating several nearby lines at once.

Focus of this adjustment

The goal is more precise control of the selected path, rather than mainly pushing the skin upward. Follow the application instructions alongside this case-specific recommendation.

05 / INITIAL SIZE COMPARISON

Compare size 2 with size 3

These are recommendations for this individual case. The technician uses the two trials to choose a starting cut size and placement.

DAY 1 / SIZE 2

Apply size 2 according to the instructions. Check placement before judging whether the size suits you. Record the requested application and response videos.

DAY 2 / SIZE 3

On the following day, repeat with size 3. Use the same recording format so the technician can compare the two trials.

What the technician compares

1. Where the fold actually forms.
2. One smooth path versus several competing lines.
3. Balance between the inner and outer corners.
4. Similarity to the existing crease on the right eyelid.
5. How the fold changes after tape removal.

View the application and removal instructions

06 / NEXT STEPS

Review the response, then adjust

The starting plan focuses on the left eyelid. The existing crease on the right remains the reference while the first two trials are reviewed.

If you also want to raise the right crease

The original videos do not show the intended higher position on the right or its response to tracing at that height. More demonstration is needed before recommending a change there.

Your next steps in this case

1. Practice controlled tracing with the single-point end.
2. Apply size 2 carefully and record the requested videos.
3. Repeat with size 3 on the next day, in the same format.
4. Send both trials for a comparison of cut size and placement.
5. Review a higher right-side crease only after showing that goal.

What we can say from this case

Faint lines, low existing creases and a temporary fold response provide starting observations. The next decision comes from the actual size-trial videos, not from a promised result date.

About this public sample

The client name and order number have been removed. Photos and case observations are retained. This detailed report illustrates individual guidance; the report and ongoing support included depend on the service selected.

Contact Optifold